

## PRODUCT SPECIFICATION

**Product:** SAVETT WOUND CLEANSER,  
REFILL  
**Reference code:** 3227



### 1. General description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantee sterility. The sachets are packed in a refill box.

### 2. Intended use

Savett Wound Cleanser is intended for mechanical cleansing of minor wounds when clean tap water is not available.

### 3. Technical description

#### *Product dimensions*

Sachet	(l x w) mm:	60 x 58 ±2
Tissue	(l x w) mm:	120 x 110 ±2
Box	(l x w x h) mm:	60 x 60 x 151

#### *Material*

Tissue:	White non-woven tissue
Foil:	Aluminum foil with print
Liquid:	Physiological saline solution

### 4. Sterilization

Sterilization method:	Irradiation
-----------------------	-------------

### 5. Composition

#### *Consumer package*

<u>Product/Material</u>	<u>Pcs</u>
Savett Wound Cleanser sachets	40

#### *Transport package*

Number of consumer packages in transport package:	10
Dimensions (l x w x h) mm:	330 x 135 x 160
Weight (kg):	1,85



## 6. Labeling

### *On product (sachet)*

Cederroth address, telephone- and  
faxnumber, web-address, mail-address  
Product name, Reference code  
Lot-number, Expiry date  
Symbols: CE +NB#, STERILE, DO NOT  
REUSE

### *On consumer package*

Cederroth address, telephone- and  
faxnumber, web-address, mail-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code  
Symbols: CE +NB#, STERILE, LOT,  
USE BY

### *On transport package*

Cederroth address, telephone- and  
faxnumber, web-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code  
Quantity  
Symbols: LOT, USE BY

## 7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

## 8. Accessories

The following accessories can be used with the product but are not necessary.

<u>Ref #</u>	<u>Product name</u>
3250	Savett Dispenser

**Authorized by:** Teija Ålander  
*Quality Assurance Manager*

Date: 2008-02-27